STERILE ULTRASOUND NEEDLES

INSTRUCTIONS FOR USE

INDICATIONS FOR USE

The ultrasound needles are intended for use with a Bausch & Lomb Microsurgical System for the emulsification of an opacified crystalline lens.

INSTRUCTIONS FOR USE: COAXIAL PHACOEMULSIFICATION

- A. Refer to FIG. 1 for sleeve adjustment. Expose about 1.0 mm of the titanium needle (A) beyond the end of the silicone sleeve (B) for most cases. This may be increased to approximately 1.5 mm for dense nuclei (see FIG. 1). For MICS needs BL3318(A)S, assemble sleeve as shown in FIG. 2. BSS wetting of the sleeve prior to assembly is advised to assist assembling the infusion sleeve to the needle.
- B. See chart for recommended size of knife to use for the incision to achieve minimum wound leakage.

Needle	Thread Color	Incision Size
DP81XX Standard	Gold	3.0 mm or greater
DP82XX MicroFlow [®] 1	Green	2.5 mm or greater
DP83XX MicroFlow [®] Plus ¹	Blue	2.85 mm or greater
DP84XX Thin Tip	Silver	2.5 mm or greater
DP87XX MicroFlow [®] 2.2 ¹	Teal	2.2 mm or greater
DP88XX Vortex 19 ga.	Silver	2.5 mm or greater
DP89XX Vortex 20 ga.	Silver	2.4 mm or greater
BL3318(A)S MICS ²	Purple	1.8 mm or greater

- C. When inserting the tip into the eye, the following technique will avoid introducing air into the anterior chamber. Enter the eye with the footpedal in position 1.
- D. During entry, the infusion sleeve portals may get hung up in the small incision. Simply rotate the handpiece back and forth, alternating pronation and supination of the wrist, in order to release the sleeve as the tip enters the eye.
- E. The infusion bottle height can be set according to physician preference. Bausch & Lomb recommends to keep the drip chamber approximately 30 inches (76.2 cm) above the patient's eye level.

INSTRUCTIONS FOR USE: BIMANUAL PHACOEMULSIFICATION

See chart below for recommended size of knife to use for the incision to achieve minimum wound leakage.

Needle	Thread Color	Incision Size
DP81XX Standard	Gold	3.0 mm or greater
DP82XX MicroFlow [®] 1	Green	2.5 mm or greater
DP83XX MicroFlow [®] Plus ¹	Blue	2.85 mm or greater
DP84XX Thin Tip	Silver	2.5 mm or greater
DP87XX MicroFlow [®] 2.2 ¹	Teal	2.2 mm or greater

- A. Attach the modified infusion sleeve, with distal tip removed, onto the phaco needle.
- B. To achieve minimum wound leakage, use a 1.4–1.6mm keratome to prepare phaco incision.
- C. The vacuum, flow rate, phaco power and bottle height could be set according to physician preference.
- Sterile ultrasound needles are for single use only, remove and discard needle from handpiece after surgery.

During any ultrasonic procedure, metal particles may result from inadvertent touching of the ultrasonic tip with a second instrument. Another potential source of metal particles resulting from any ultrasonic handpiece may be the result of ultrasonic energy causing micro abrasion of the ultrasonic tip.

^{*/} are trademarks of Bausch & Lomb Incorporated or its affiliates.
© Bausch & Lomb Incorporated. All rights reserved.

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